


Body Analyser TANITA BC-418MA

SPECIFICATIONS

The body analyser is designed to determine the weight and percentage of body fat, an indicator of the healthy range of body fat, etc., see the overview of measurable parameters / output information.

OVERVIEW OF MEASURABLE PARAMETERS/ OUTPUT INFORMATION

- Total weight [kg]
- BMI, muscle tissue weight [kg]
- Total body water [kg]
- BMR Basal metabolism [kJ, kcal and BMR Indicator]
- Body fat weight [kg] and body fat percentage
- Visceral fat indicator
- Healthy range of body fat
- Weight without fat [kg]
- Muscle mass [kg]
- Bone mass weight [kg]
- Bone mineral weight [kg]
- Segmental evaluation:
 - Body fat [kg, %]
 - Fat free mass (FFM) [kg]
 - Muscle mass [kg]



TYPES OF SAMPLES SUITABLE FOR ANALYSIS / MEASUREMENT CONDITIONS

Man / woman aged 18-99 years.